GEBRAUCHSANWEISUNG
Coaxial I/A System Disposable - DE

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INSTRUCTIONS FOR USE
Coaxial I/A System Disposable - EN

INTENDED USE

The Coaxial (mono-manual) I/A Disposable Instrument is used for a combination of delivering solutions to the eye and extracting solutions and residuals of tissue from the eye during cataract operations in the anterior chamber. The instrument is connected to the aspiration / irrigation system by two standard Luer-cones.

The Bimanual I/A Disposable Instruments are intended to be used only by appropriately specialist surgeon.

IMPORTANT USER INFORMATION

Insufficient flow conditions can cause an imbalance of the fluids. It is strongly recommended to observe the pressure conditions and, if necessary, to adjust them. It is very important to prevent any collapse of the anterior chamber. Bürki inno med AG disclaims any liability for any inappropriate and incorrect handling.

CHECK BEFORE USE

Following inspection are necessary before use:

- Visual inspection of obstruction, deformation and no cracks
- Filling of irrigation and aspiration system with NaCl 0.9% solution to ensure the functionality

WARNINGS/ PRECAUTIONS (Regard packaging symbols, too.)

- The Coaxial (mono-manual) I/A Disposable Instrument is designed for single use; any re-sterilisation or re-use is strongly forbidden and can lead to dangerous post-operative infections.
- Don't use the instrument after expiry date or if the packaging is damaged.
- Prior the use, inspect the tip of the Coaxial (mono-manual) I/A System. Don't use the product, if the tip is damaged.
- The instrument is only to be used by qualified professionals.
- The manufacturer shall not be liable for any injury or damage suffered by a patient due the use of the product.
- The Bimanual I/A System are packed in a Blister and a Tyvek-Bag. The Blister represents the sterile barrier. The Tyvek-Bag is an additional protection which is intended to facilitate the handling in operating room.

TECHNICAL SPECIFICATION

Parameter	Specification
Life span	< 60min
Connectors	Standard Luer-cones
Shelf life	3 years in original packaging
Sterilisation type	Gamma sterilisation
State of delivery	Sterile
Quality Management System of Bürki inno med AG	EN ISO 13485

DISPOSAL

After usage, the Bimanual I/A Disposable Instruments must be handled as a contaminated item.

The devices must be disposed according to local regulations. For information, contact your local environmental or public office and appropriate waste disposal companies. In this regard, waste is to be recycled or disposed of:

- without danger to human health
- without using procedures or methods harmful to the environment, particularly to the water, air, soil, flora and fauna
- without creating noise or smell

DEVICE LABEL AND SYMBOLS

	Single use only	REF	Article number		Do not resterilize
	Do not use after expiry date	STERILE R	Sterilisation by irradiation		Flow direction row
LOT	Batch / Lot Number		Consult instructions for use		Manufacturer
EC REP	Authorised representative Bürki inno med GmbH Im Schaffner 47/1, 69123 Heidelberg-Pfaffengrund, Germany		Do not use if the packaging has been damaged	CE 0297	ID number of notified body
	Single sterile barrier system	Manufacturer:	Bürki inno med AG · Industriestrasse 67 · CH-9443 Widnau · Switzerland Tel. +41 (0)71 744 88 66 · www.buerki.		

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Bürki inno med AG rejette toute responsabilité pour les applications inappropriées.