

## Vitreotomy Disposable Forceps & Scissors



### Product Features:

- Single use / no risk of cross contamination
- Double packed, Gamma - sterilisation
- Lightweight
- Manual pressure for handle ~ 220g or ~150g
- Plastic handle / Stainless steel tip
- 23G / 25G
- End gripping micro forceps tips
- Serrated micro forceps tip
- Pack size of 10
- Units per box: 50

### Forceps in details



### Available in January 2017:

- |          |  |
|----------|--|
| VDF 2311 | Standard disposable micro forceps; serrated jaws 23G |
| VDF 2312 | Standard disposable micro forceps; end gripping 23G  |

### Available in March 2017:

- |          |  |
|----------|--|
| VDF 2511 | Standard disposable micro forceps; serrated jaws 25G |
| VDF 2512 | Standard disposable micro forceps; end gripping 25G  |

Version 08/2016