

Coaxial Disposable Irrigation Aspiration System 20G, 21G, 22G




Product Features:

- Coaxial System Disposable 20G & 22G, 21G optional
- Lightweight irrigation / aspiration Handpiece
- Prolonged sleeve (plus 1.5mm)
- all G: 45° tip smooth abrasive (rough) or polished
- 20G, 21G: 65° tip smooth abrasive (rough) or polished
- Compatible with your existing Phaco or I/A System
- Swiss Design / Swiss Made


Details:




prolonged tip



45° tip



65° tip



./.

CX2045-P

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 45°, POLISHED

CX2065-P

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 65°, POLISHED

CX2145-P

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 45°, POLISHED

CX2165-P

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 65°, POLISHED

CX2245-P

Coaxial Sytem Disposable 22G
Asp. port 0.35mm, angled 45°, POLISHED

CX2045-R

Coaxial System Disposable 20G
Asp. Port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)

CX2065-R

Coaxial System Disposable 20G
Asp. port 0.35mm, angled 65°, SMOOTH
ABRASIVE (ROUGH)

CX2145-R

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)

CX2165-R

Coaxial System Disposable 21G
Asp. port 0.35mm, angled 65°, SMOOTH
ABRASIVE (ROUGH)

CX2245-R

Coaxial Sytem Disposable 22G
Asp. port 0.35mm, angled 45°, SMOOTH
ABRASIVE (ROUGH)