

Bimanual Reusable Family Aspiration 20G, 21G, 22G, 23G


Product Features:

- Lightweight, colour coded aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made




Aspiration

20G




BA 108/T special model
Aspiration handle smooth
abrasive (rough) tip,
(0.80mm), Asp. Port 0.35mm


21G




BA 100/T
Aspiration Handle polished
tip, (0.80mm); Asp. Port
0.30mm




BA 102/T
Aspiration Handle polished tip,
(0.80mm); Asp. Port 0.40mm



BA 104/T
Aspiration Handle polished
tip, (0.80mm); Asp. Port
0.50mm



BA 107/T special model
Aspiration Handle smooth
abrasive (rough) tip,(0.80mm);
Asp. Port 0.40mm




BA 105/T
Aspiration Handle smooth
abrasive (rough)
tip,(0.80mm); Asp. Port
0.30mm




BA 118/T
Aspiration Handle polished
tip, (0.70mm); Asp. Port
0.35mm

BA 119/T special model
Aspiration Handle smooth
abrasive (rough) tip,
(0.70mm); Asp. Port 0.35mm

22G




BA 113/T
Aspiration Handle polished
tip, (0.65mm); Asp. Port
0.30mm



BA 114/T
Aspiration Handle smooth
abrasive (rough) tip,(0.65mm);
Asp. Port 0.30mm

23G



Version 05/2022

About Bürki

Bürki plans, develops and produces
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
info@buerki.eu
www.buerki.eu