

Reusable Bimanual Irrigation Aspiration System 20G, 21G



Product Features:

- Lightweight, colour coded irrigation and Aspiration Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made

Most often used Bimanual I/A Systems

Aspiration


21G




BA 100/T
Aspiration Handle polished tip, 21G (0.80mm); Asp. Port 0.30mm



BA 102/T
Aspiration Handle polished tip, 21G (0.80mm); Asp. Port 0.40mm



BA 104/T
Aspiration Handle polished tip, 21G (0.80mm); Asp. Port 0.50mm



BA 105/T
Aspiration Handle smooth abrasive (rough) tip 21G (0.80mm); Asp. Port 0.30mm




BA 107/T Special model
Aspiration Handle smooth abrasive (rough) tip 21G, (0.80mm); Asp. Port 0.40mm

Irrigation



BI 120/T
Irrigation Handle polished tip, 21G (0.80mm), dual side ports 0.50mm




BI 125/T
Irrigation Handle smooth abrasive (rough) tip, 21G (0.80mm), dual side ports 0.50mm



BI 140/T
Irrigation Handle polished tip, 21G (0.80mm), open end 45°

20G




BA 108/T special model
Aspiration Handle smooth
abrasive (rough) tip 20G,
(0.90mm); Asp. Port
0.35mm




BI 126/T
Irrigation Handle open end tip
45°, 20G (0.90mm), High flow


Oval: based on 21G




BAO 100/T
Aspiration Handle oval
cannula, 0.55mm x
0.95mm, mate finish tip;
Asp. Port 0.30mm



BIO 120/T
Irrigation Handle oval cannula,
0.55mm x 0.95mm, mate finish
tip, Irrig. Port 0.55mm



BAO 102/T special model
Aspiration Handle oval
cannula, 0.70mm x
1.10mm, smooth abrasive
(rough); Asp. Port 0.35mm



BIO 122/T
Irrigation Handle oval cannula,
0.55mm x 0.95mm, mate finish
tip, Irrig. Port 0.45mm



BAO 105/T
Aspiration Handle oval
cannula, 0.55mm x
0.95mm, smooth abrasive
(rough) tip;
Asp. Port 0.30mm



BIO 127/T
Irrigation Handle oval cannula,
1.10mm x 0.70mm, mate finish,
open end tip 30°, High flow