

Bimanual Reusable Family Irrigation 20G, 21G, 22G, 23G

Product Features:

- Lightweight, colour coded irrigation Handpieces
- Compatible with your existing Phaco or I/A System
- Smooth abrasive (rough) tip available for capsule polishing
- Swiss Design / Swiss Made



Irrigation

20G



BI 126/T
Irrigation Handle open end tip
45° (0.90mm), High flow

21G



BI 120/T
Irrigation Handle polished tip
(0.80mm), dual side ports
0.50mm



BI 125/T
Irrigation Handle smooth
abrasive (rough) tip,
(0.80mm), dual side ports
0.50mm



BI 140/T
Irrigation Handle polished tip
open end 45°

22G



BI 128/T special model
Irrigation Handle polished tip,
dual side ports 0.45mm



BI 129/T special model
Irrigation Handle polished tip,
open end 45°

23G



BI 123/T
Irrigation Handle polished tip
(0.65mm), dual side ports
0.35mm



BI 124/T
Irrigation Handle polished tip
(0.65mm), open end 45°



BI 127/T special model (MLI)
Irrigation Handle polished tip
(0.65mm)

Version 05/2022

About Bürki
Bürki plans, develops and produces
microsurgical instruments in SWISS QUALITY
EN ISO 13485

Bürki inno med AG
Industriestrasse 67
9443 Widnau
Switzerland

Tel +41 (0)71 744 88 66
info@buerki.eu
www.buerki.eu